

# ONLY ONE

**FL<sub>OW</sub>SAFE<sup>®</sup> II**  
DISPOSABLE CPAP SYSTEM **EZ**

## CPAP System with Integrated Nebulizer.

A major leap in product innovation,  
it takes Emergency Care to a whole new level.

- CPAP and Nebulizer System in One
- Less Parts -- Easy Set-up
- **Only One** O<sub>2</sub> Source



**MercuryMedical<sup>®</sup>**  
Your Need . . . Our Innovation<sup>®</sup>

[www.mercurymed.com](http://www.mercurymed.com)

# FLOW-SAFE<sup>®</sup> II

DISPOSABLE CPAP SYSTEM

## EZ

When lives hang in the balance, you need a disposable CPAP system that lives up to the challenge.

The Flow-Safe<sup>®</sup> II EZ Disposable CPAP System from Mercury Medical is the only system that has so many advantages, it gives you an advantage in life-saving emergencies.

In short, the Flow-Safe II EZ represents nothing short of a revolution in emergency care.

**SAFE.**  
**SURE.**  
**SUPERIOR.**  
**CPAP.**



Clinicians can deliver aerosol & CPAP therapy with just one oxygen source.

CPAP system uses 50% less oxygen,  
Increase flow and activate the nebulizer.



**ADVANTAGE:**  
INTEGRATED NEBULIZER

The Flow-Safe II EZ offers the capability of an in-line nebulizer.

- Uses only one oxygen supply source.
- Easy Set-Up, less parts.

**ADVANTAGE:**

CPAP SYSTEM CONSUMES 50% LESS OXYGEN

The Flow-Safe II EZ conserves oxygen while maintaining high FiO<sub>2</sub> delivery. Saves resources and represents a major advantage on long transports. Increasing flow may be necessary when activating the nebulizer.



**ADVANTAGE:**

BUILT-IN MANOMETER AND PRESSURE RELIEF VALVE

Only the Flow-Safe II EZ system features a built-in manometer for verified pressure readings. No assembly of separate apparatus. And the pressure relief valve automatically adjusts to avoid excess pressure.

**ADVANTAGE:**

ADVANCED MASK DESIGN

The Flow-Safe II EZ mask is lighter, easier to handle and is designed to form a better anatomical seal. The elastic head harness is easy-to-place with velcro straps that easily adjust for patient comfort.

- Includes Quick Disconnect Clips
- Straight Rotating Port
- Soft Forehead Padding easily adjusts reducing pressure on the nose.



Clinicians can *now* deliver aerosol & CPAP therapy with just one integrated oxygen source.

CPAP system uses 50% less oxygen,

Increase flow and activate the nebulizer.



CO<sub>2</sub> Sampling Port  
Fits In-Line  
(Supplied Separately)

*Ordering Information*

Part Number	Description	Quantity
10-57318	Flow-Safe II EZ, Large Adult Deluxe Mask with Ports and EZ Flow Max Nebulizer	5/Box
10-57319	Flow-Safe II EZ, Small Adult Deluxe Mask with Ports and EZ Flow Max Nebulizer	5/Box
10-57320	Flow-Safe II EZ, Child Deluxe Mask with Ports and EZ Flow Max Nebulizer	5/Box
10-57321	Flow-Safe II EZ, without Mask with EZ Flow Max Nebulizer	5/Box

**FASTER. SAFER. MORE COMFORTABLE.  
MORE CONVENIENT. MORE EFFICIENT.**